

Maximized Exposure.
Optimized Protection.

SurgiSleeve[™] Wound Protector



Maximized Wound Exposure

Circumferential elastic retraction maximizes working area for better exposure and visualization.





Product Details



The proximal blue ring is easy to roll down to adjust to variable thicknesses of the abdominal wall.

3x stronger material offers optimized wound protection¹

The flexible distal grey ring is inserted through the incision in seconds and easily removed at the end of the procedure.

ALWAYS REFER TO PRODUCT INSTRUCTIONS FOR USE FOR COMPLETE INDICATIONS, CONTRAINDICATIONS, WARNINGS, PRECAUTIONS AND OPERATING INSTRUCTIONS.

Sizes and Reorder Codes

Reorder Codes	Description	Quantity per Box
WPSM256	Wound Protector Small Incision Size 2.5 - 6cm	5
WPMD509	Wound Protector Medium Incision Size 5 - 9cm	5
WPLG914	Wound Protector Large Incision Size 9 - 14cm	5

Puncture resistance comparison of mean strengths between SurgiSleeve™ wound protector and Alexis™* wound protector as part of a February 9, 2012 Puncture Resistance Material Evaluation (Report # 2151-002). Statistically significant results p=0.004.



IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

COVIDIEN, COVIDIEN with logo, Covidien logo and "positive results for life" are U.S. and internationally registered trademarks of Covidien AG.

TM * Trademark of its respective owner. Other brands are trademarks of a Covidien company. © 2013 Covidien. - M130534 - 05/2013



IO CONTACT US, PLEASE VISIT WWW.COVIDIEN.COM AND SELECT "CONTACT US